

PIEZO ULTRASONIC SCALERS

25% OFF

D600 LED

- Smart touching technique system.
- Clear control panel.
- Ultra-sensitive.
- Waterproof design.
- Aluminium Shell Handpiece
- Portable design with water reservoir.

Reg. \$1399.99

NET
\$1049⁹⁹

3 Modes:

1. Supragingival Scaling
2. Subgingival Scaling
3. Endodontic Irrigation

ONE YEAR WARRANTY



S6 LED

- Clear control panel display,
- Ultrasensitive touch control.
- Waterproof design.
- Aluminium Shell Handpiece

Reg. \$1599.99

NET
\$1199⁹⁹

3 Modes:

1. Supragingival Scaling
2. Subgingival Scaling
3. Endodontic Irrigation

ONE YEAR WARRANTY



DTE D5 SCALER (SATELEC² TYPE)

Scaling, Perio, Endo modes

DTE-D5 SCALER
HANDPIECE AND TIPS
ARE COMPATIBLE WITH
SATELEC² TYPE SCALERS

DTE-D5 LED Piezo Ultrasonic Scaler
complete with 6 Tips
ED1, GD1, GD2, GD3, GD4, PD1
and Male Quick Connector

Reg. \$899.99

NET
\$674⁹⁹

ONE YEAR WARRANTY
(except Scaler Tips)



DETACHABLE HANDPIECE

- Digitally controlled, automatic frequency tracking
- Small, light weight and easy operation
- Autoclavable Titanium handpiece with excellent piezo-ceramic vibrations
 - Frequency: 28kHz ±3kHz



DTE D7 (SATELEC² TYPE) PORTABLE SCALER

Scaling, Perio, Endo modes

**NO FIXED
INSTALLATION**

Reg. \$1199.99

NET
\$899⁹⁹

DTE D7 PORTABLE SCALER
HANDPIECE AND TIPS
ARE COMPATIBLE WITH
SATELEC² TYPE SCALERS

DTE-D7 LED Portable Scaler
complete with 6 Tips

GD1, GD4, GD5, GD6, PD1, PD3 and PD4



1 EMS is a registered trademark of Electro Medical System S.A. Switzerland
2 SATELEC is a registered trademark of Satelec/Acteon Group France

SCALING

- G1T** (WOODPECKER, EMS¹, Type)
GD1T (DTE, SATELEC², NSK³ Type)
for supragingival calculus and
plaque, neck of teeth and hard
calculus on the tooth point
each **\$34⁹⁹**
- G2T** (WOODPECKER, EMS¹, Type)
GD2T (DTE, SATELEC², NSK³ Type)
for supragingival heavy calculus
each **\$34⁹⁹**
- G3T** (WOODPECKER, EMS¹, Type)
GD3T (DTE, SATELEC², NSK³ Type)
for subgingival, supragingival and
interdental calculus and plaque
each **\$34⁹⁹**
- G4T** (WOODPECKER, EMS¹, Type)
GD4T (DTE, SATELEC², NSK³ Type)
for supragingival calculus and plaque
each **\$34⁹⁹**
- G5T** (WOODPECKER, EMS¹, Type)
GD5T (DTE, SATELEC², NSK³ Type)
for supragingival and neck calculus
and plaque
each **\$39⁹⁹**
- G6T** (WOODPECKER, EMS¹, Type)
GD6T (DTE, SATELEC², NSK³ Type)
for supragingival heavy calculus
each **\$39⁹⁹**

WOODPECKER SCALER TIPS STAINLESS STEEL TITANIUM COATED

Non-Sticky, Longer Lasting

For WOODPECKER, EMS¹, DTE, SATELEC², NSK³ Type Scalers

PERIO

- P1T** (WOODPECKER, EMS¹, Type)
PD1T (DTE, SATELEC², NSK³ Type)
for subgingival and supragingival
calculus
each **\$39⁹⁹**
- P3T** (WOODPECKER, EMS¹, Type)
PD3T (DTE, SATELEC², NSK³ Type)
for fine scaling and deep perio calculus
each **\$39⁹⁹**
- P4T** (WOODPECKER, EMS¹, Type)
PD4T (DTE, SATELEC², NSK³ Type)
for subgingival and interdental calculus
each **\$39⁹⁹**

C&B

- G7T** (WOODPECKER, EMS¹, Type)
GD7T (DTE, SATELEC², NSK³ Type)
used to remove dental crown
each **\$44⁹⁹**

ENDO

- E1T** (WOODPECKER, EMS¹, Type)
ED1T (DTE, SATELEC², NSK³ Type)
(120° Angle Holder) for
NITI K-Files, usually used for
anterior root canal cleaning,
sterilizing and debris removal
each **\$44⁹⁹**
- E2T** (WOODPECKER, EMS¹, Type)
ED2T (DTE, SATELEC², NSK³ Type)
(95° Angle Holder) for
NITI K-Files, usually used for
anterior root canal cleaning,
sterilizing and debris removal
each **\$44⁹⁹**
- E3T** (WOODPECKER, EMS¹, Type)
ED3T (DTE, SATELEC², NSK³ Type)
After root canal cleaning, using
the mechanism and heat energy of
ultrasonic vibration to intertense the
filling lateral condensation, to make
sure the top of the scaler tip touch the
lateral condensation
each **\$44⁹⁹**

Buy 4 get 1 FREE
(same item)

1 EMS is a registered trademark of Electro Medical System S.A. Switzerland
2 SATELEC is a registered trademark of Satelec/Acteon Group France
3 NSK is a registered trademark of Nakanishi Inc. Japan