



safco Brand Products

Low Prices. Quality Service.

Excel Impression Material

Excel VPS Imp. Material
Buy 4, Get 1 FREE!

4x 50ml Cart.	Regular	Fast
Monophase	023-1020524	023-1020559 (1/mo)
Light Body	023-1020502	023-1020533 (3/mo)
Heavy Body	023-1020518	023-1020547 (6/mo)

\$44.99 \$35.99NET

023-1020309 Empty Storable Canister **\$10.99**



Nice Bite Bite Registration

2 x 50ml Cart.

023-3410024 Superfast Set (4/mo) **\$34.99**
023-3410013 Fast Set (8/mo)



Excel Alginate

Similar to Kromopan

023-1020304 Excel Alginate Chroma 11lb Pouch **\$12.99**



10 + 2 FREE! OR 16+4 FREE!

Surgical Aspirators

032-1022719 White 1/8" Hole, 6.5" Length (25) (10/mo)
032-1022733 Green 1/4" Hole, 6.5" Length (25) (14/mo)
032-1022715 Blue 1/16" Hole, 7.75" Length (25) (14/mo)

\$9.99



6 + 2 FREE! OR 20+10 FREE!

Kleer Mirrors

032-1023850 Cone Socket #4 (10/mo)
032-1023857 Cone Socket #5 (15/mo)
032-1023864 Simple Stem #4 (2/mo)
032-1028870 Simple Stem #5 (3/mo)

\$29.99 12/box



3 + 1 FREE!

T-Mixers

Acrylic C&B T-Mixer mixing tips, 25/pkg



025-1022463 T-Mixer (9/mo) 10:1 Mixing Tip (C&B)

\$29.99

4 + 1 FREE!



025-1020588 Small Yellow (30/mo) (Light Body) (15/mo)

025-1020597 Large Teal (Heavy Body / Bite Reg)

\$49.99

Intraoral Syringe Tips

Buy 2, Get 1 FREE!

085 Yellow for yellow mixing tips
086 Clear for Green or Pink mixing tips

\$34.99 \$23.33NET



Saliva Ejector

- Side vents with narrow, recessed center suction hole.
- Not made with natural rubber latex.

10 + 2 FREE!

40+10 FREE!



\$4.99 \$3.99NET

032-1027600 Clear With Blue Tip (100) (300/mo)
032-1027656 White With White Tip (100) (100/mo)

Chill Refrigerant Spray 10 oz

(15/mo)

4 + 1 FREE!



Cold stimulus spray for the rapid detection of pulp vitality. Provides a low temperature of -26.2°C.

Wintergreen scent.

Clear 4" dispensing straw included.

Similar to Endo-Ice.

016-1028550 **\$34.99 \$27.99NET**

Nice Bite Dual Arch Tray



Dual Arch Bite Trays, 48/pk

025-1025350 Anterior
025-1025354 Posterior

\$39.99

40% more volume

safco Brand Products

Low Prices. Quality Service.



Tartar & Stain Remover

Ready to use

019-0815661 1 Gal. (10-21/mo)

\$29.99

3 + 1 FREE!



Be-Safe

4 + 1 FREE!

Ultrasonic Cleaner Tablets

019-1020192 64/box (15/mo)

\$49.99



Enzymatic Evacuation Cleaner

032-1020074 1.89L (NON FOAMING) 1:32

\$24.99

3 + 1 FREE!



Instrument Ties

Autoclavable up to 275°F.

Colours: Blue, Green, Purple, Red, Yellow

6/pk \$29.99

032-4051701 - 09



3 + 1 FREE!

Sonic-Clean

Ultrasonic Cleaner

\$699.95

019-1020172

Safco Sonic-Kleen features durable design and construction and delivers 43kHz of consistent, dependable power year after year.

Includes Basket

Tank capacity:	0.85 Gallon
Use capacity:	3.2L
Overall Dimensions:	10.3" x 6.5" x 8.3"
Tank Dimensions:	9.4" x 13.7" x 4"
Output:	95 watts / 0.8 amps
Operation Frequency:	43 kHz



Kleen & Safe Headrest Covers

4 + 1 FREE!



011-0364976 11.25" x 10" (250/box) \$29.99

011-0364972 9.5" x 11" (250/box) \$24.99

Tray Sleeve

4 + 2 FREE!




032-1022675 10.5" x 14" (500/box) \$41.99

STD TRAY "B" \$27.99NET

Mixing Pads

4 + 1 FREE!

100 Sheets



025-0486240 1" x 1" (#0, X-Small) \$2.99

025-0486245 3" x 3" (#1, Small) \$3.99

025-0486250 3" x 6" (#3, Large) \$5.99

025-0486255 6" x 6" (#4, X-Large) \$9.99

Mouth Props

4 + 1 FREE!

2/pkg	Rubber	Silicone
Small Child	032-0498645	032-0498655
Child	032-0498646	032-0498656
Adult	032-0498647	032-0498657

\$24.99



Multi-Purpose Dispensing Needle Tips

3 + 1 FREE!

20+10 FREE!



100/pk \$12.99

(Note: Tip colour may vary)

010-1026103	18GA. (LINER) PINK HUB
010-1026107	20GA. (FLOW) BLACK HUB
010-1026110	20GA. (FLOW) YELLOW HUB
010-1026113	22GA. GREY HUB
010-1026115	22GA. (ETCH) BLUE HUB

Matrix Bands

3 + 1 FREE!

10+5 FREE!

Thickness 0.002 Thickness 0.0015

004-1063312 #1 004-1063340 #1

004-1063351 #2 004-1063363 #2

004-1063324 #3 004-1063378 #3

004-1063335 #13 004-1063389 #13

\$7.99 36/pk

CONTOUR MATRIX BAND

004-1025519 Ultrathin \$47.99 50/pk
004-1025515 Regular

4 + 1 FREE!





safco Brand Products

Low Prices. Quality Service.

Labstone Buff

- Cured compressive strength: 8,500 psi
- Initial set time: 10-12 minutes
- Mixing Ratio: 30ml/100g
- Setting expansion: 0.12%

102-0415 Buff, 33 lb. carton

\$54.99



Lead Free Aprons

- Cured compressive strength: 8,500 psi • Initial set time: 10-12 minutes
- Mixing Ratio: 30ml/100g • Setting expansion: 0.12%

039-1020386 Adult, with Collar - Beige **\$299.99**

039-1020382 Adult, with Collar - Blue


039-1020389 Child, with Collar - Blue White, Dog Paw Pattern

039-1020394 Thyroid Protector - Beige **\$79.99**

039-1020390 Thyroid Protector - Blue

Panoramic Adult Aprons w/o Collar - Beige/Blue **\$399.99**

Panoramic Child Aprons w/o Collar - Beige/Blue **\$349.99**



SafeWrap™ Basics

X-Ray Syringe Sleeves



- Convenient opening for easy insertion.
- Precut slit for air/water syringe tip.
- 2.5" x 10"

102-2670 Clear, box of 500

\$19.99

4 + 1 FREE!

20+10 FREE!

X-Pose Ultra D58 Single Film



4 + 1 FREE!

9 + 3 FREE!

039-1027559 D-58 #2 Film (1 x 150/box) **\$59.99**

Universal Film Mounts Tear Apart

Holds #0, #1, or #2 films vertically or horizontally. Roll of 1,000 mounts in dispenser box

\$44.99



4 + 1 FREE!

Safe-Wrap Basics Barrier film

1200 Sheets

019-1022778 Clear **\$22.99**

019-1022773 Blue

3 + 1 FREE!



Safco Polyback Bibs

2ply paper + 1 poly



500pcs

\$34.99

\$27.99NET

3 + 1 FREE!

Colours: Light Blue, Lilac, Green, White.

Paper Cups Polycoated



2 + 1 FREE!

4 oz. (1,000/case)

011-0206604 Bubbles

011-0206606 NEW Green

\$49.99 \$33.33NET

5 oz. (1,000/case)

011-0206605 Bubbles

011-0206607 NEW Lavender

\$51.99 \$34.66NET

Ritter Tray Cover

1000/case, White

011-0801499 8.5"x12.25"

3 + 1 FREE!

\$34.99 \$26.24NET



LARR-TEST

24 Hours 100/box

BIOLOGICAL INDICATOR FOR STEAM AUTOCLAVE

019-0082233



1@ **\$209.99**

3@ **\$199.99**